

Virus Filtration in Continuous Bioprocessing – Considerations for Filter Design Space and Validation Strategies

Morven McAlister¹, Aernout Martens², Nigel Jackson³ and Stijn de Backer⁴

¹25 Harbor Park Drive, Port Washington, USA; ²Nijverheidsweg, Medemblik, The Netherlands; ³5 Harbourgate Business Park, Portsmouth, UK; ⁴2 Reugelstraat, Hoegaarden, Belgium

INTRODUCTION

Continuous bioprocessing of biopharmaceuticals can offer many benefits such as more agile manufacturing with a reduced footprint. Chromatography and virus inactivation are often implemented as continuous unit operations, and viral clearance is validated by similar strategies to a batch process. However, most end-users have leaned towards batch-mode implementation of virus filters, even in a continuous process, by repeatedly storing product in a tank, prior to passing through a virus filter, and consequently batch-mode virus filter validation.

However, to fully realize the benefits of continuous bioprocessing, it is desirable to fully integrate virus filtration into continuous bioprocessing applications. Pall have applied the principles of Quality by Design (QbD) to assess the different design spaces for batch vs. continuous bioprocessing and identify potential risks. By analyzing both supplier and end-user data, we have developed technical solutions and are proposing three different strategies that can be used to overcome some of the complexities associated with validating virus filtration in a fully integrated continuous bioprocess.

INTEGRATION OF VIRUS FILTRATION (VF) INTO CONTINUOUS BIOPROCESSING – ASSESSMENT OF DESIGN SPACE

Figure 1




Figure 2






Figure 3

Compilation of Pall and >50 end-user virus validation results supports robust performance of Pegasus™ Prime virus filters in continuous bioprocessing design space



Through regulatory engagement and prior knowledge assessments, low flux, extended filtration times and process interruptions are considered CPPs for integrating virus filtration into continuous bioprocessing applications. However, multiple studies have demonstrated robust performance of Pegasus Prime virus filters at the extremes of these CPPs.

This novel design space presents additional complexities for validation of virus filtration into a fully integrated continuous bioprocess. Therefore, Pall have developed solutions to overcome some of these perceived barriers, and can offer strategies to address different technical challenges.



PROPOSED VIRUS FILTER VALIDATION STRATEGIES FOR INTEGRATED CONTINUOUS BIOPROCESSING APPLICATIONS

Figure 4
PP7 titer is maintained in different fluids for up to 22 days




Figure 5
MVM titers in cell culture media are not impacted at any temperature tested for 48 hours




Figure 6
X-MuLV titers in cell culture media are negatively impacted by time and temperature




Figure 7




Figure 8






Figure 9



- ▶ Pre-testing is recommended to verify specific virus titer in product over time prior to performing virus filtration studies
- ▶ Pall have proposed three different validation strategies for integration of virus filtration into continuous bioprocessing
 - **One Spike – One Study** – can be used if test virus maintains titer over test duration. Most simplistic virus filter validation strategy.
 - **Spike Replenishment** – spike and process fluid are changed at pre-determined frequencies. Requires additional sampling and adds complexity to determine overall virus clearance.
 - **Inline Spiking** – consider for situations when viruses are removed by the prefilter and/or the mAb reagggregates after prefiltration. Requires careful consideration of flow dynamics to ensure adequate mixing of the virus under low flow conditions.



CONCLUSIONS

- ▶ Technical complexities associated with complete integration of virus filtration into continuous bioprocessing can be perceived as a barrier to implementation
- ▶ Pall have demonstrated that by following QbD principles to assess CPPs, alternative validation strategies are required
- ▶ Next generation filters such as the Pegasus Prime virus filter, show robust performance in the unique design space required for continuous bioprocessing
- ▶ Determining a suitable virus filter validation strategy requires careful consideration of virus and product stability, but Pall have developed three methods to address different scenarios
- ▶ Alternative validation strategies were based on regulatory input for integrated continuous bioprocessing
- ▶ Pall is available for further technical assistance with all virus filtration applications

International Offices

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.